Lord's

LordsMed

Urine Analyzer
LM 1023

10 And 14 Parameters Strip Options Inbuilt Printer

Better supported for outside testing



One Touch Operation

Small in size Touch screen

Portable | **High-efficient** | **Smart**

Lord's

LordsMed

Features



Easy to operate
One-click Operation, click "start" to complete automatically

Rapid
14 Items in 60s

Tiny and Portable
178 mm 87 mm 46 mm, 360g

Intelligent Calibration
Automatic Color Compensation

Large Battery Capacity
4-hour Continuous Working

Scenarios


Laboratory


Clinics


Outdoor Rescue


Hospitalization&E.R


Family doctor


Ambulance

Lord's

LordsMed

Function



Operations



Step.1 Sample Loading

Step.2 Reaction

Step.3 Determination

Parameters of The Device

Device	Urine Analyzer
Detection principle	Dual Wavelength Reflection Method
Sample type	Urine
Size	178 mm*87 mm*46 mm
Weight	360g
Test Item	LEU URO mAlb PRO BIL GLU ASC SG KET NIT CRE pH BLO CA